

## BagSeal®

Thermal sealing unit

Ref 261 000



Manufactured under:



BagSeal® allows a wide and tight sealing of all kind of bags (with/without filter). Perfect closure is guaranteed.

### TECHNICAL SPECIFICATIONS

- All stainless steel except handles, welding blades (anodized aluminum), bar joint (chromated galvanized steel)
- Large sealing: 250 x 5 mm
- Adjustable power heating
- LED light code at the end of the sealing (welding process duration indicator: welding + cooling time)
  
- For bags 80/100/400/2000 mL
- Pressure: 22 Newton/cm<sup>2</sup>
  
- Voltage and frequency: 100-127 V~ / 220-240 V~ 50-60 Hz
- Max power: 1600 W
- International protection marking: 40
  
- Dimensions (w x d x h): 35.9 x 27.6 x 32.2 cm, weight: 10.8 kg
- Box (w x d x h): 49 x 40 x 49 cm, gross weight: 13.6 kg
  
- Warranty 1 year
  
- Manufactured under CE, RoHS, WEEE
- Designed and made in France

### OPTIONS / ACCESSORIES

- Ref 261 500 Maintenance kit for BagSeal®: 2 sealing threads + 2 fuses + 2 Teflon strips

#### PARIS

Phone: +33 1 34 62 62 61  
[info@interscience.com](mailto:info@interscience.com)

#### FRANKFURT

Phone: +49 611 7238 7770  
[sales.germany@interscience.com](mailto:sales.germany@interscience.com)

#### BOSTON

Phone: +1 781 937 0007  
[sales.usa@intersciencelab.com](mailto:sales.usa@intersciencelab.com)

#### SHANGHAI

电话: +86 178 2123 6642  
[sales.china@interscience.com](mailto:sales.china@interscience.com)

#### SINGAPORE

Phone: +65 6977 7232  
[sales.asia@interscience.com](mailto:sales.asia@interscience.com)

#### TOKYO

Phone: +81 3 6712 9715  
[sales.japan@interscience.com](mailto:sales.japan@interscience.com)